

Product Name: Recombinant Human MRPL34 Protein
Catalog #: PEH10810



Summary

| | |
|------------------------------|--|
| Name | RM34_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9BQ48 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 6 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | KARGNEYQPS NIKRKNKHGW VRRLLSTPAGV QVILRRMLKG RKLSLH |

image



Background

Alternative Names 39S ribosomal protein L34, mitochondrial (L34mt) (MRP-L34) (Mitochondrial large ribosomal subunit protein bL34m)

Note

For Research Use Only , Not for Diagnostic Use