


Product Name: Recombinant Human GADD45A Protein
Catalog #: PEH10815



Summary

Name	GA45A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P24522
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKTERMDKVG DALEEVLSKA LSQRTITVG V YEAAKLLNVD PDNVVLCCLA ADEDDDRDVA LQIHFTLIQA FCCENDINIL RVSNPGR LAE LLLLET DAGP AASEGAEQPP DLHCVLVTNP HSSQWKDPAL SQLICFCRES RYMDQWVPVI

image



Background

Alternative Names	Growth arrest and DNA damage-inducible protein GADD45 alpha (DNA damage-inducible transcript 1 protein) (DDIT-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use