


Product Name: Recombinant Human MSX1 Protein
Catalog #: PEH10816



Summary

Name	MSX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P28360
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FGKPAGGGAG QAPSAAAATA AAMGADEEGA KPKVSPSLLP FSVEALMADH RKPGAKESAL APSEGVQAAG GSAQPLGVPP GSLGAPDAPS SPRPLGHFSV GLLKLPEDA LVKAESPEKP E RTPWMQSPR

image



Background

Alternative Names Homeobox protein MSX-1 (Homeobox protein Hox-7) (Msh homeobox 1-like protein)

Note

For Research Use Only , Not for Diagnostic Use