


Product Name: Recombinant Human TAL1 Protein
Catalog #: PEH10820



Summary

Name	TAL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17542
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSPYEMEITD GPHTKVVRRI FTNSRERWRQ QNVNGAFAEL RKLIPHTPPD KKLSKNEILR LAMKYINFLA KLLNDQEEEG TQRAKTGKDP VVGAGGGGGG GGGAPPDDL LQDVLSPNSS CGSSLDGAAS

image



Background

Alternative Names T-cell acute lymphocytic leukemia protein 1 (TAL-1) (Class A basic helix-loop-helix protein 17) (bHLHa17) (Stem cell protein) (T-cell leukemia/lymphoma protein 5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TAL1 Protein
Catalog #: PEH10820

