


Product Name: Recombinant Human ETV6 Protein
Catalog #: PEH10823



Summary

Name	ETV6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P41212
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HTQPEVILHQ NHEEDNCVQR TPRPSVDNVH HNPPTIELLH RSRSPITTNH RSPDPEQRP LRSPLDNMIR RLSPAERAQG PRPHQENNHQ ESYPLSVSPM ENNHCPASSE SHPKPSSPRQ ESTRVIQLMP SPIMHPLILN PRHSVDFKQS RLSEDGLHRE GKPINLSHRE

image



Background

Alternative Names	Transcription factor ETV6 (ETS translocation variant 6) (ETS-related protein Tel1) (Tel)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ETV6 Protein
Catalog #: PEH10823

