


Product Name: Recombinant Human DYRK1A Protein
Catalog #: PEH10826



Summary

Name	DYR1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13627
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFPAPLGWSG TEAPTQVTVE THPVQETTFH VAPQQNALHH HHGNSSHHHH HHHHHHHHHG QQALGNRTRP RVYNSPTNSS STQDSMEVGH SHHSMTSLSS STTSSTSSS STGNQGNQAY QNRPVAANTL DFGQNGAMDV NLTVYSNPRQ ETGIAGHPTY

image



Background

Alternative Names	Dual specificity tyrosine-phosphorylation-regulated kinase 1A (EC 2.7.12.1) (Dual specificity YAK1-related kinase) (HP86) (Protein kinase minibrain homolog) (MNBH) (hMNB)
--------------------------	--

Note

Product Name: Recombinant Human DYRK1A Protein
Catalog #: PEH10826



For Research Use Only , Not for Diagnostic Use