

Product Name: Recombinant Human MBD4 Protein
Catalog #: PEH10836



Summary

Name	MBD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95243
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DDVNFRKVRK PKGKVTILKG IPIKTKKGC RKSCSGFVQS DSKRESVCNK ADAESEPVAQ KSQLDRTVCI SDAGACGETL SVTSEENSLV

image



Background

Alternative Names	Methyl-CpG-binding domain protein 4 (EC 3.2.2.-) (Methyl-CpG-binding endonuclease 1) (Methyl-CpG-binding protein MBD4) (Mismatch-specific DNA N-glycosylase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use