


**Product Name: Recombinant Human CDK8 Protein**  
**Catalog #: PEH10844**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | CDK8_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | P49336   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 14 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | PTSDVFAGCQ IPYPKREFLT EEPDDKGDK KNQQQQGNN HTNGTGHPGN<br>QDSSHTQGPP LKKVRVPPT TSGGLIMTS DYQRSNPAA YPNPGPSTSQ<br>PQSSMGYSAT  |

## image



## Background

|                          |  |
|--------------------------|--|
| <b>Alternative Names</b> | Cyclin-dependent kinase 8 (EC 2.7.11.22) (EC 2.7.11.23) (Cell division protein kinase 8) (Mediator complex subunit CDK8) (Mediator of RNA polymerase II transcription subunit CDK8) (Protein kinase K35) |
|--------------------------|--|

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human CDK8 Protein**  
**Catalog #: PEH10844**

---

