


Product Name: Recombinant Human SP4 Protein
Catalog #: PEH10845



Summary

Name	SP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02446
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EAAAAAAMAT EGGKTSEPEN NNKKPKTSGS QDSQPSPLAL LAATCSKIGT PGENQATGQQ QIIIDPSQGL VQLQNQPQQL ELVTTQLAGN AWQLVASTPP ASKENNVSQP ASSSSSSSSS NNGSASPTKT KSGNSSTPGQ FQVIQVQNPS GSVQYQVIPQ LQTVGEQQIQ INPTSSSSLQ DLQGQIQLIS AGNNQAILTA ANRTASGNIL AQNLANQTVP VQIRPGVSIP

image



Background

Alternative Names Transcription factor Sp4 (SPR-1)

Note

Product Name: Recombinant Human SP4 Protein
Catalog #: PEH10845



For Research Use Only , Not for Diagnostic Use