



Product Name: Recombinant Human AFF1 Protein
Catalog #: PEH10846



Summary

Name	AFF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P51825
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERRNQEAHQE KEAFPEKIPL FGEPYKTAKG DELSSRIQNM LGNYEEVKEF LSTKSHTHRL DASENRLGKP KYPLIPDKGS SIPSSSFHTS VHHQSIHTPA SGPLSVGNIS HNPCKMAQPRTEPMPSLHAKSCGPPDSQHLT QDRLGQEGFG SSHHKKGDDR

image



Background

Alternative Names	AF4/FMR2 family member 1 (ALL1-fused gene from chromosome 4 protein) (Protein AF-4) (Protein FEL) (Proto-oncogene AF4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AFF1 Protein
Catalog #: PEH10846

