


Product Name: Recombinant Human REV1 Protein
Catalog #: PEH10848



Summary

Name	REV1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBZ9
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLSFNPVCRP EDPLGPSNI AKQLNNRVNH IVKKIETENE VKVNGMNSWN EEDENND FSF VDLEQTSPGR KQNGIPHPRG STAIFNGHTP SSNGALKTQD CLVPMVNSVA SRLSPAFSQE EDKAEKSSTD FRDCTLQQLQ QSTRNTDALR NPHRTNSFSL SPLHSNTKIN GAHHSTVQGP SSTKSTSSVS TFSKAAPSVP SKPSDCNFIS

image



Background

Alternative Names	DNA repair protein REV1 (EC 2.7.7.-) (Alpha integrin-binding protein 80) (AIBP80) (Rev1-like terminal deoxycytidyl transferase)
--------------------------	---

Note

Product Name: Recombinant Human REV1 Protein
Catalog #: PEH10848



For Research Use Only , Not for Diagnostic Use