


Product Name: Recombinant Human TCEA1 Protein
Catalog #: PEH10850



Summary

Name	TCEA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P23193
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEPAITSQNS PEAREESTSS GNVSNRKDET NARDTYVSSF PRAPSTSDSV RLKCREMLAA ALRTGDDYIA IGADEEELGS QIEEAIYQEI RNTDMKYKNR VRSRISNLKD AKNPNLRKNV LCGNIPDLF ARMTAEEMAS DELKEMRKNL TKEAIREHQM AKTGGTQTDL

image



Background

Alternative Names	Transcription elongation factor A protein 1 (Transcription elongation factor S-II protein 1) (Transcription elongation factor TFIIS.o)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TCEA1 Protein
Catalog #: PEH10850

