


Product Name: Recombinant Human HAND2 Protein
Catalog #: PEH10851



Summary

Name	HAND2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P61296
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGVPPGAGPP GLGGPRPVKR RGTANRKERR RTQSINSAFA ELRECIPNVP ADTKLSKIKT LRLATSYIAY LMDLLAKDDQ NGEAEAFKAE IKKTDVKEEK RKKELNEILK

image



Background

Alternative Names	Heart- and neural crest derivatives-expressed protein 2 (Class A basic helix-loop-helix protein 26) (bHLHa26) (Deciduum, heart, autonomic nervous system and neural crest derivatives-expressed protein 2) (dHAND)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HAND2 Protein
Catalog #: PEH10851

