


Product Name: Recombinant Human HSPA9 Protein
Catalog #: PEH10854



Summary

Name	GRP75_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P38646
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDRRKKERVE AVNMAEGIIH DTETKMEEFK DQLPADECNK LKEEISKMRE LLARKDSETG ENIRQAASSL QQASLKLFEF AYKKMASERE GSGSSGTGEQ

image



Background

Alternative Names Stress-70 protein, mitochondrial (75 kDa glucose-regulated protein) (GRP-75) (Heat shock 70 kDa protein 9) (Mortalin) (MOT) (Peptide-binding protein 74) (PBP74)

Note

For Research Use Only , Not for Diagnostic Use