


Product Name: Recombinant Human TGIF1 Protein
Catalog #: PEH10855



Summary

Name	TGIF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15583
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVESVMGIKN FMPALEETPF HSCTAGPNPT LGRPLSPKPS SPGSVLARPS VICHTTVTAL KDVPFSLCQS VGVGQNTDIQ QIAAKNFTDT SLMYPEDTCK

image



Background

Alternative Names Homeobox protein TGIF1 (5'-TG-3'-interacting factor 1)

Note

For Research Use Only , Not for Diagnostic Use