


Product Name: Recombinant Human VSX1 Protein
Catalog #: PEH10859



Summary

Name	VSX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZR4
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APSRPPPALG RQKRSDSVST SDEDSQSEDR NDLKASPTLG KRKKRRHRTV FTAHQLEELE KAFSEAHYPD VYAREMLAVK

image



Background

Alternative Names Visual system homeobox 1 (Homeodomain protein RINX) (Retinal inner nuclear layer homeobox protein) (Transcription factor VSX1)

Note

For Research Use Only , Not for Diagnostic Use