


Product Name: Recombinant Human SIN3A Protein
Catalog #: PEH10862



Summary

Name	SIN3A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96ST3
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HQGINFKQND TKVLRSKSLL NEIESIYDER QEQATEENAG VPVGPHLSLA YEDKQILED AALIIHHVKR QTGIQKEDKY KIKQIMHHFI PDLLFAQRGD LSDVEEEEE EMDVDEATGA VKKHNGVGGS PPKSKLLFSN TAAQKLRGMD EVYNLFYVNN NWYIFMRLHQ ILCLRLRIC SQAERQIEEE NREREWEREV LGIKRDKSDS PAIQLRLKEP MDVDVEDYYP AFLDMVRSLL DGNIDSSQYE DSLREMFTIH AYIAFTMDKL IQSIVRQLQH IVSDEICVQV TDLYLAENNN GATGGQLNTQ NSRSLLESTY QRKAEQLMSD ENCFKLMFIQ SQGQVQLTIE LLDTEENSD DPVEAERWSD YVERYMNSDT

image



Background

Alternative Names	Paired amphipathic helix protein Sin3a (Histone deacetylase complex subunit)
--------------------------	--

Product Name: Recombinant Human SIN3A Protein
Catalog #: PEH10862



Sin3a) (Transcriptional corepressor Sin3a)

Note

For Research Use Only , Not for Diagnostic Use