


## Summary

<b>Name</b>	KRT34_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & N-His
<b>Accession #</b>	O76011
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	54 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	QLTCFSITCS STMSYSCCLP SLGCRTSCSS RPCVPPSCHG YTLPGACNIP ANVSNCNWFC EGSFNGSEKE TMQFLNDRLA SYLEKVRQLE RDNAELEKLI QERSQQQEPL LCPSYQSYFK TIEELQQKIL CAKAENARLV VNIDNAKLAS DDFRSKYQTE QSLRLLVESD INSIRRILDE LTLCKSDLES QVESLREELI CLKKNHEEEV NTLRSQLGDR LNVEVDTAPT VDLNQVLNET RSQYEALVEI NRREVEQWFA TQTEELNKQV VSSSEQLQSC QAEIHELRRRT VNALEIELQA QHNLRDLEN TLTESEAHYS SQLSQVQSLI TNVESQLAEI

## image



## Background

**Alternative Names** Keratin, type I cuticular Ha4 (Hair keratin, type I Ha4) (Keratin-34) (K34)

**Product Name: Recombinant Human KRT34 Protein**  
**Catalog #: PEH10863**

---



**Note**

For Research Use Only , Not for Diagnostic Use