

Product Name: Recombinant Human KRT34 Protein
Catalog #: PEH10864



Summary

Name	KRT34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O76011
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPCVPPSCHG YTLPGACNIP ANVSNCNWFC EGSFNGSEKE TMQFLNDRLA SYLEKVRQLE RDNAELEKLI QERSQQQEPL LCPSYQSYFK TIEELQQKIL CAKAENARLV VNIDNAKLAS DDFRSKYQTE QSLRLLVESD INSIRRLILDE LTLCKSDLES QVESLREELI CLKKNHEEEV NTLRSQLGDR LNVEVDTAPT VDLNQVLNET RSQYEALVEI NRREVEQWFA TQTEELNKQV VSSSEQLQSC QAEIIELRRT VNALEIELQA QHNLRDLEN TLTESEAHYS SQLSQVQSLI TNVESQLAEI RCDLERQNQE YQVLLDVRAR LECEINTYRS

image



Background

Alternative Names Keratin, type I cuticular Ha4 (Hair keratin, type I Ha4) (Keratin-34) (K34)

Product Name: Recombinant Human KRT34 Protein
Catalog #: PEH10864



Note

For Research Use Only , Not for Diagnostic Use