


Product Name: Recombinant Human ONECUT2 Protein
Catalog #: PEH10865



Summary

Name	ONEC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95948
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HPHHHPHHHH HHHHQRLSGN VSGSFTLMRD ERGLPAMNNL YSPYKEMPGM SQLSPLAAT PLGNLGLLH NAQQLPNYG PPGHDKMLSP NFDAHHTAML TRGEQHLSRG LGTPPAAMMS HLNGLHHPGH TQSHGPVLAP

image



Background

Alternative Names One cut domain family member 2 (Hepatocyte nuclear factor 6-beta) (HNF-6-beta) (One cut homeobox 2) (Transcription factor ONECUT-2) (OC-2)

Note

For Research Use Only , Not for Diagnostic Use