


Product Name: Recombinant Human ZNF169 Protein
Catalog #: PEH10871



Summary

Name	ZN169_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14929
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EALMAFRDVA VAFTQKEWKL LSSAQRTLYR EVMLENYSHL VSLGIAFSKP KLIEQLEQGD EPWRENEHL LDLCPEPRTE FQPSFPHLVA FSSQLLRQY ALSGHPTQIF PSSSAGGDFQ LEAPRCSEK GESGETEGPD SSLRKRPRI SRTFFSPHQG DPVEWVEGMR EGGTDLRLAQ RMSLGGSDTM LKGADTSESG

image



Background

Alternative Names Zinc finger protein 169

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF169 Protein
Catalog #: PEH10871

