


Product Name: Recombinant Human SAFB Protein
Catalog #: PEH10872



Summary

Name	SAFB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15424
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSSCGRNFVW SGLSSTTRAT DLKNLFSKYG KVVGAKVVTN ARSPGARCYG FVTMSTAEAA TKCINHLHKT ELHGKMISVE KAKNEPVGKK TSDKRDSGK

image



Background

Alternative Names Scaffold attachment factor B1 (SAF-B) (SAF-B1) (HSP27 estrogen response element-TATA box-binding protein) (HSP27 ERE-TATA-binding protein)

Note

For Research Use Only , Not for Diagnostic Use