


Product Name: Recombinant Human TMEM119 Protein
Catalog #: PEH10876



Summary

Name	TM119_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q4V9L6
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQKASAYYPS SFPKKKYVDQ SDRAGGPRAF SEVPDRAPDS RPEEALDSSR QLQADILAAT QNLKSPTRAA LGGGDGARMV EGRGAEETEEK GSQEGDQEVQ GHGVPVETPE AQEPCSGVL

image



Background

Alternative Names Transmembrane protein 119 (Osteoblast induction factor) (OBIF)

Note

For Research Use Only , Not for Diagnostic Use