


Product Name: Recombinant Human ZNF470 Protein
Catalog #: PEH10878



Summary

Name	ZN470_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ECI4
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GIKLCKAMSL GSVTFDVAI DFSQDEWEWL NLAQRSLYKK VMLENYRNLV SVGLCISKPD VISLLEQEKD PWVIKGGMNR GLCPDLECVW VTKSLSLNQD IYEEKLPPAI IMERLKSVDL ECSTLGKNWK CEDLFEELV NQKTHFRQET ITHIDTLIEK RDHSNKSSTV FHLNTLSYIK QIFPMEERIF NFHTDKKSLK THSVVKKHKQ DRGEKLLKC NDCEKIFSKI STLTLHQRIH TGEKPYECIE

image



Background

Alternative Names Zinc finger protein 470 (Chondrogenesis zinc finger protein 1) (CZF-1)

Note

Product Name: Recombinant Human ZNF470 Protein
Catalog #: PEH10878



For Research Use Only , Not for Diagnostic Use