


Product Name: Recombinant Human NHLRC1 Protein
Catalog #: PEH10880



Summary

Name	NHLC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6VVB1
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGWGTLVNPT GLALCPKTGR VVVVHDGRRR VKIFDSGGGC AHQFGEKGDA AQDIRYPVDV TITNDCHVVV TDAGDRSIKV FDFGQIKLV IGGQFSLPWG VETTPQNGIV VTDAEAGSLH LLDVDFAEV LRRTERLQAH LCNPRGVAVS WLTGAIIVLE HPLALGTGVC STRVKVFSS MQLVGQVDTF GLSLYFPSKI TASAVTFDHQ GNVIVADTSG PAILCLGKPE EFPVPKPMVT HGLSHPVALT FTKENLLVL DTASHSIKVY

image



Background

Alternative Names E3 ubiquitin-protein ligase NHLRC1 (EC 2.3.2.27) (Malin) (NHL repeat-containing protein 1) (RING-type E3 ubiquitin transferase NHLRC1)

Product Name: Recombinant Human NHLRC1 Protein
Catalog #: PEH10880



Note

For Research Use Only , Not for Diagnostic Use