


Product Name: Recombinant Human DDX46 Protein
Catalog #: PEH10882



Summary

Name	DDX46_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L014
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KAEGKIIKKS SGFSGKGFKF DETEQALANE RKKLQKAALG LQSDDEDAA VDIDEQIESM FNSKKRVKDM AAPGTSSVPA PTAGNAEKLE IAKRLALRIN AQKNLGIESQ DVMQQATNAI LRGGTILAPT VSAKTIAEQL AEKINAKLNY VPLEKQEEER QDGGQNESFK RYEEELEIND FPQTARWKVT SKEALQRIS YSEAAITIRG TYFPPGKEPK EGERKIYLAI ESANELAVQK AKAEITRLIK EELIRLQNSY

image



Background

Alternative Names	Probable ATP-dependent RNA helicase DDX46 (EC 3.6.4.13) (DEAD box protein 46) (PRP5 homolog)
--------------------------	--

Note

Product Name: Recombinant Human DDX46 Protein
Catalog #: PEH10882



For Research Use Only , Not for Diagnostic Use