


Product Name: Recombinant Human ZDHHC17 Protein
Catalog #: PEH10886



Summary

Name	ZDH17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IUH5
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MADGPDEYDT EAGCVPLLHP EEIKPQSHYN HGYGEPLGRK THIDDYSTWD IVKATQYGIY ERCRELVEAG YDVRQPDKEN VTLLHWAAIN NRIDLKYYI SKGAIVDQLG GDLNSTPLHW ATRQGHLSMV VQLMKYGADP SLIDGEGCSC IHLAAQFGHT SIVAYLIAKG QDVDMMDQNG MTPLMWAAYR THSVDPTRLL LTFNVSVNLG DKYHKNTALH WAVLAGNTTV ISLLEAGAN VDAQNIKGES ALDLAKQRKN VWMINHLQEA RQAKGYDNPS FLRKLKADKE

image



Background

Alternative Names	Palmitoyltransferase ZDHHC17 (EC 2.3.1.225) (Acyltransferase ZDHHC17) (EC 2.3.1.-) (DHHC domain-containing cysteine-rich protein 17) (DHHC17) (Huntingtin yeast partner H) (Huntingtin-interacting protein 14) (HIP-14)
--------------------------	---

Product Name: Recombinant Human ZDHHC17 Protein
Catalog #: PEH10886



(Huntingtin-interacting protein 3) (HIP-3) (Huntingtin-interacting protein H)
(Putative MAPK-activating protein PM11) (Putative NF-kappa-B-activating
protein 205) (Zinc finger DHHC domain-containing protein 17)

Note

For Research Use Only , Not for Diagnostic Use