


**Product Name: Recombinant Human TMTC1 Protein**  
**Catalog #: PEH10888**



## Summary

<b>Name</b>	TMTC1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8IUR5
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	48 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LFRSGVQTLP HNAKVHYNIA NFLKDQGRNK EAIYHYRTAL KLYPRHASAL NNLGLTLTRDT AEAKMYYQRA LQLHPQHNRA LFNLGNNLLKS QEKKEEAITL LKDSIKYGPE FADAYSSLAS LLAEQERFKE AEEIYQTGIK NCPDSSDLHN NYGVFLVDTG LPEKAVAHYQ QAIKLSPSHH VAMVNLGRLY RSLGENSMAE EWYKRALQVA HKAEILSPLG ALYYNTGRYE EALQIQEAA ALQPSQRELR LALAQVLAVM GQTKEAEKMT NHIVSEETGC LECYRLLSAI YSKQENHDKA LDAIDKALQL KPKDPKVISE LFFTKGNQLR EQNLLDKAFE SYRVAVQLNP DQAQAWMNMG GIQHIKGKYV SARAYYERL QLVPDSKLLK ENLAKLDRLE KRLQEVREKD

## image



## Background

**Product Name: Recombinant Human TMTC1 Protein**  
**Catalog #: PEH10888**



---

**Alternative Names**

Protein O-mannosyl-transferase TMTC1 (EC 2.4.1.109) (Transmembrane and TPR repeat-containing protein 1)

**Note**

For Research Use Only , Not for Diagnostic Use