


**Product Name: Recombinant Human CSGALNACT2 Protein**  
**Catalog #: PEH10889**



## Summary

|                              |                                                                                                                                                                                                                                                                                      |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | CGAT2_HUMAN                                                                                                                                                                                                                                                                          |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                     |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                |
| <b>Accession #</b>           | Q8N6G5                                                                                                                                                                                                                                                                               |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                               |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                |
| <b>Molecular Mass</b>        | 30 kDa                                                                                                                                                                                                                                                                               |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                            |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                       |
| <b>Sequence</b>              | NASLPGVVGE NYGKEYYQAL LQEQEEHYQT RATSLKRQIA QLKQELQEMS<br>EKMRS LQERR NVGANGIGYQ SNKEQAPSDL LEFLHSQIDK AEVSIGAKLP<br>SEYGVIPFES FTLMKVFQLE MGLTRHPEEK PVRKDKRDEL VEVIEAGLEV<br>INNPDEDDEQ EDEEGPLGEK LIFNENDFVE GYYRTERDKG TQYELFFKKA<br>DLTEYRHVTL FRPFGPLMKV KSEMIDITRS IINIIVPLAE |

## image



## Background

|                          |                                                                                                                                                                      |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Chondroitin sulfate N-acetylgalactosaminyltransferase 2 (EC 2.4.1.174)<br>(Chondroitin beta-1,4-N-acetylgalactosaminyltransferase 2)<br>(Beta4GalNAcT-2) (GalNAcT-2) |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|



**Product Name: Recombinant Human CSGALNACT2 Protein**  
**Catalog #: PEH10889**

---

**Note**

For Research Use Only , Not for Diagnostic Use