


Product Name: Recombinant Human RNF128 Protein
Catalog #: PEH10895



Summary

Name	RN128_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TEB7
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRKQRQLKAD AKKAIGRLQL RTLKQGDKEI GPDGDSCAVC IELYKPNDLV RILTCNHIFH KTCVDPWLL E HRTCPMCKCD ILKALGIEVD VEDGSVSLQV PVSNEISNSA SSHEEDNRSE TASSGYASVQ GTDEPPLEEH VQSTNESLQL VNHEANSVAV DVIPHVDNPT FEEDETPNQE

image



Background

Alternative Names	E3 ubiquitin-protein ligase RNF128 (EC 2.3.2.27) (Gene related to anergy in lymphocytes protein) (GRAIL) (RING finger protein 128) (RING-type E3 ubiquitin transferase RNF128)
--------------------------	--

Note

Product Name: Recombinant Human RNF128 Protein
Catalog #: PEH10895



For Research Use Only , Not for Diagnostic Use