


Product Name: Recombinant Human SAMD8 Protein
Catalog #: PEH10903



Summary

Name	SAMD8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96LT4
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RWTTKHVAVW LKDEGFFEYV DILCNKHLRD GITLLTLEY DLRSPPLEIK VLGDIKRLML SVRKLQKIHI DVLEEMGYNS DSPMGSMTPF ISALQSTDWL CNGELSHDCD GPITDLNSDQ YQYMNGKNKH SVRRLDPEYW

image



Background

Alternative Names	Sphingomyelin synthase-related protein 1 (SMSr) (EC 2.7.8.-) (Ceramide phosphoethanolamine synthase) (CPE synthase) (Sterile alpha motif domain-containing protein 8) (SAM domain-containing protein 8)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SAMD8 Protein
Catalog #: PEH10903

