



Product Name: Recombinant Human IWS1 Protein
Catalog #: PEH10904



Summary

Name	IWS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96ST2
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKNEEKDLFG SDSESGNEEE NLIADIFGES GDEEEEFETG FNQEDLEEEK GETQVKEAED SDSDDNIKRG KHMDFLSDFE MMLQRKKSMS GKRRRNRDGG TFISDADDVV SAMIVKMNEA AEEDRQLNNQ KKPALKKLTL LPAVVMHLKK QDLKETFIDS GVMSAIKEWL SPLPDRSLPA LKIREELKI LQELPSVSQE TLKHSGIGRA VMYLYKHPKE SRSNKDMAGK LINEWSRPIF GLTSNYKGMT REEREQRDLE

image



Background

Alternative Names Protein IWS1 homolog (IWS1-like protein)

Product Name: Recombinant Human IWS1 Protein
Catalog #: PEH10904



Note

For Research Use Only , Not for Diagnostic Use