


Product Name: Recombinant Human ZDHHC5 Protein
Catalog #: PEH10905



Summary

| | |
|------------------------------|--|
| Name | ZDHC5_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9C0B5 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 38 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SSAKLSRGDS LKEPTSIAES SRHPSYRSEP SLEPESFRSP TFGKSFHFDP LSSGSRSSSL KSAQGTGFEL GQLQSIRSEG TTSTSYKSLA NQTRNGSLSY DSSLTPSDSP DFESVQAGPE PDPPLGYTSP FLSARLAQQR EAERHPRLVP TGPTHREPSP VRYDNLSRHI VASLQEREKL LRQSPPLPGR EEEPGLGDSG IQSTPGSGHA PRTSSSSDDS KRSPLGKTPL GRPAVPRFGK PDGLRGRGVG SPEPGTAPY LGRSMSYSSQ KAQPGVSETE EVALQPLLTP KDEVQLKTTY SKSNGQPKSL GSASPGPGQP PLSSPTRGGV KKVSGVGGTT |

image



Background

Alternative Names Palmitoyltransferase ZDHHC5 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 5) (DHHC-5) (Zinc finger protein 375)

Product Name: Recombinant Human ZDHHC5 Protein
Catalog #: PEH10905



Note

For Research Use Only , Not for Diagnostic Use