


**Product Name: Recombinant Human ST6GALNAC1 Protein**  
**Catalog #: PEH10911**



## Summary

<b>Name</b>	SIA7A_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9NSC7
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	14 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	KPSRHQRTEK IKERSLQSLA KPKSQAPTRA RRTTIYAEPV PENNALNTQT QPKAHTTGDR GKEANQAPPE EQDKVPHTAQ RAAWKSPEKE KTMVNTLSPR GQDAGMASGR

## image



## Background

<b>Alternative Names</b>	Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 1 (EC 2.4.99.3) (GalNAc alpha-2,6-sialyltransferase I) (ST6GalNAc I) (ST6GalNAcI) (Sialyltransferase 7A) (SIAT7-A)
--------------------------	--

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human ST6GALNAC1 Protein**  
**Catalog #: PEH10911**

---

