


Product Name: Recombinant Human RBM28 Protein
Catalog #: PEH10912



Summary

Name	RBM28_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NW13
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEAAKLQTTK VKKPTGTRNL YLAREGLIRA GTKAAEGVSA ADMAKRERFE LLKHQKLDQ NIFVSRLC LHNLPKAVDD KQLRKL LLSA TSGEKG VRIK ECRVMRDLKG VHG NMGQSL GYAFAEFQEH EHALKALRLI NNNPEIFGPL KRPIVEFSLE DRRKLKMKEL RIQRSLQKMR SKPATGEPQK GQPEPAKDQQ QKAAQHHTEE QSKVPPEQKR KAGSTSWTGF QTKAEVEQVE LPDGKKRRKV LALPSHRGPK IRLRDKGKVK PVHPKPKPQ INQWKQEKQQ LSSEQVSRKK AKGNKTETRF NQLVEQYKQK LLGPSKGAPL

image



Background

Alternative Names RNA-binding protein 28 (RNA-binding motif protein 28)

Product Name: Recombinant Human RBM28 Protein
Catalog #: PEH10912



Note

For Research Use Only , Not for Diagnostic Use