



Product Name: Recombinant Human MYEF2 Protein
Catalog #: PEH10913



Summary

Name	MYEF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9P2K5
Host	E.coli
Species	Human
Molecular Mass	77 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPAEPPGEPR REPHPAEAEK QQPQHSSSSN GVKMENDESA KEEKSDLKEK STGSKKANRF HPYSKDKNSG AGEKKGPNRN RVFISNIPYD MKWQAIKDLM REKVG EVTIV ELFKDAEGKS RCGVVEFKD EEVKKALET MNKYDLSGRP LNIKEDPDGE NARRALQRTG GSFPGGHVDP MGSGLMNLPP SILNNPNIPP EVISNLQAGR LGSTIFVANL DFKVGWKKLK EVFSIAGTVK RADIKEDKDG KSRGMGTVTF EQAIEAVQAI SMFNGQFLFD RPMHVKMDDK SVPHEEYRSH

image



Background

Alternative Names	Myelin expression factor 2 (MEF-2) (MyEF-2) (MST156)
--------------------------	--

Product Name: Recombinant Human MYEF2 Protein
Catalog #: PEH10913



Note

For Research Use Only , Not for Diagnostic Use