


Product Name: Recombinant Human STX18 Protein
Catalog #: PEH10914



Summary

Name	STX18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2W9
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDGSRDELFR RSPRPGDFS SRAREVISHI GKLRDFLLEH RKDYINAYSH TMSEYGRMTD TERDQIDQDA QIFMRTCSEA IQQLRTEAHK EIHSQQVKEH RTAVLDFIED YLKRVCKLYS EQRAIRVKRV VDKKRLSKLE PEPNTKTRES TSSEKVSQSP SKDSEENPAT EERPEKILAE TQPELGTWGD GKGEDELSPE EIQMFEQENQ RLIGEMNSLF DEVRQIEGRV VEISRLQEIF TEKVLQQEAE IDSIHQLVVG ATENIKEGNE

image



Background

Alternative Names	Syntaxin-18 (Cell growth-inhibiting gene 9 protein)
--------------------------	---

Product Name: Recombinant Human STX18 Protein
Catalog #: PEH10914



Note

For Research Use Only , Not for Diagnostic Use