


Product Name: Recombinant Human XPR1 Protein
Catalog #: PEH10915



Summary

Name	XPR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBH6
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITPEWRKQYI QYEAFKDMLY SAQDQAPSVE VTDEDTVKRY FAKFEEKFFQ TCEKELAKIN TFYSEKLAEA QRRFATLQNE LQSSLDAQKE STGVTTLRQR RKPVFHLSHE ERVQHRNIKD LKLAFSEFYI SLILLQNYQN LNFTGFRKIL KKHDKILETS RGADWRVAHV EVAPFYTCKK INQLISETEA VVTNELEDGD RQKAMKRLRV PPLGAAQPAP

image



Background

Alternative Names Xenotropic and polytropic retrovirus receptor 1 (Protein SYG1 homolog) (Xenotropic and polytropic murine leukemia virus receptor X3) (X-receptor)

Note

Product Name: Recombinant Human XPR1 Protein
Catalog #: PEH10915



For Research Use Only , Not for Diagnostic Use