


Summary

Name	ABCD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBJ2
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEVTELAGYT ARVYNMFVWF DEVKRGYKR TAVIQEESH SKNGAKVELP LSDTLAIK GK VIDVDHGIC ENVPIITPAG EVVASRLNFK VEEGMHLLIT GPNGCGKSSL FRILSGLWPV YEGVLYKPPP QHMFYIPQRP YMSLGSLRDQ VIYPDSVDDM HDKGYTDQDL ERILHNVHLY HIVQREGGWD AVMDWKDVLS GGEKQRMGMA RMFYHKPKYA LLDICTSAVS IDVEGKIFQA AKGAGISLLS ITHRPSLWKY HTHLLQFDGE GGWRFEQLDT AIRLTLSEEK QKLESQLAGI PKMQQLNEL CKILGEDSVL KTIKNEDETS

image



Background

Alternative Names	ATP-binding cassette sub-family D member 2 (Adrenoleukodystrophy-like 1) (Adrenoleukodystrophy-related protein) (hALDR)
--------------------------	-------------------------------------------------------------------------------------------------------------------------

Product Name: Recombinant Human ABCD2 Protein
Catalog #: PEH10916



Note

For Research Use Only , Not for Diagnostic Use