

Product Name: Recombinant Human PCDHGB3 Protein
Catalog #: PEH10919



Summary

Name	PCDGF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5G1
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPIRYAIP EE LDRGSLVGNL AKDLGFGVGD LPTRNLRVIA EKKFFTVSPE NGNLLVSDRI DREEICGKKS TCVLEFEMVA EKPLNFFHVT VLIQDINDNP PTFSQNITEL EISELALTGA TFALESAQDP DVGVNSLQQY YLSPDPHFSL IQKENLDGSR YPELVLKAPL DREEQPHHHL VLTAVDGGEP SRSC TTQIRV IVADANDNPP VFTQDMYRVN VAENLPAGSS VLKVM AIDMD EGINAEI IYA FINIGKEVRQ LFKLDSKTGE LTTIGELDFE ERDSYTIGVE AKDGGGHHTAY CKVQIDISDE

image



Background

Alternative Names Protocadherin gamma-B3 (PCDH-gamma-B3)

Product Name: Recombinant Human PCDHGB3 Protein
Catalog #: PEH10919



Note

For Research Use Only , Not for Diagnostic Use