


Product Name: Recombinant Human PRDX4 (N-6His)
Catalog #: PEH1310



Summary

Name	Peroxiredoxin-4/PRDX4
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Peroxiredoxin-4 is produced by our E.coli expression system and the target gene encoding Trp38-Asn271 is expressed with a 6His tag at the N-terminus.
Accession #	Q13162
Host	E.coli
Species	Human
Predicted Molecular Mass	28.9 KDa
Formulation	Supplied as a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	Peroxiredoxin-4; Antioxidant Enzyme AOE372; AOE37-2; Peroxiredoxin IV; Prx-IV; Thioredoxin Peroxidase A0372; Thioredoxin-Dependent Peroxide Reductase A0372; PRDX4
--------------------------	--

Product Name: Recombinant Human PRDX4 (N-6His)
Catalog #: PEH1310



Background

Peroxiredoxin-4 (PRDX4) is a member of the AhpC/TSA family. PRDX4 is a cytoplasmic protein and contains one thioredoxin domain. PRDX4 exists in homodimer or heterodimer with PRDX1. PRDX4 reduces hydrogen peroxide and alkyl hydroperoxides to water and alcohol with the use of reducing equivalents derived from thiol-containing donor molecules. In addition, PRDX4 is probably involved in redox regulation of the cell, regulating the activation of NF-kappa-B in the cytosol by a modulation of I-kappa-B-alpha phosphorylation.

Note

For Research Use Only , Not for Diagnostic Use.