

Product Name: Recombinant Human PFDN4 (N-6His)
Catalog #: PEH1313




Summary

Name	PFDN4/PFD4
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Prefoldin Subunit 4 is produced by our E.coli expression system and the target gene encoding Met1-Ser134 is expressed with a 6His tag at the N-terminus.
Accession #	Q9NQP4
Host	E.coli
Species	Human
Predicted Molecular Mass	17.48 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human PFDN4 (N-6His)
Catalog #: PEH1313



Alternative Names

Prefoldin Subunit 4; Protein C-1; PFDN4; PFD4

Background

Prefoldin Subunit 4 (PFDN4) is a heterohexameric chaperone protein that belongs to the prefoldin subunit beta family. The complex of PFDN4, consisting of two PFD-alpha type and four PFD-beta type subunits, forms a double beta barrel assembly with six protruding coiled-coils. PFDN4 binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly.

Note

For Research Use Only , Not for Diagnostic Use.