


Summary

Name	Protachykinin-1/TAC1
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Protachykinin-1 is produced by our E.coli expression system and the target gene encoding Glu20-Arg129 is expressed with a 6His tag at the N-terminus.
Accession #	P20366
Host	E.coli
Species	Human
Predicted Molecular Mass	15.2 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, 2M Urea, 1mM EDTA, pH 8.0.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image

Product Name: Recombinant Human TAC1 (N-6His)
Catalog #: PEH1380



Alternative Names

Protachykinin-1; PPT; TAC1

Background

Protachykinin-1(TAC1) is a secreted protein and belongs to the tachykinin family. TAC1 is encoded by the TAC1 gene. This gene encodes four products of the tachykinin peptide hormone family, substance P and neurokinin A, as well as the related peptides, neuropeptide K and neuropeptide gamma. These hormones are thought to function as neurotransmitters which interact with nerve receptors and smooth muscle cells. They are known to induce behavioral responses and function as vasodilators and secretagogues.

Note

For Research Use Only , Not for Diagnostic Use.