

Product Name: Recombinant Human SNCG
Catalog #: PEH1545




Summary

Name	SNCG/Gamma-synuclein
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Gamma-Synuclein is produced by our E.coli expression system and the target gene encoding Met1-Asp127 is expressed.
Accession #	O76070
Host	E.coli
Species	Human
Predicted Molecular Mass	13.6 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human SNCG
Catalog #: PEH1545



Alternative Names

Gamma-Synuclein; Breast Cancer-Specific Gene 1 Protein; Persyn; Synoretin; SR; SNCG; BCSG1; PERSYN; PRSN

Background

Gamma-Synuclein (SNCG) is a member of the Synuclein protein family. Gamma-Synuclein is mostly expressed in the peripheral nervous system and retina. Gamma-Synuclein plays a role in neurofilament network integrity and may be involved in modulating axonal architecture during development and in the adult. In addition, it may also function in modulating the keratin network in skin. SNCG expression in breast tumors has been as a marker for tumor progression. SNCG is also believed to be involved in the pathogenesis of neurodegenerative diseases.

Note

For Research Use Only , Not for Diagnostic Use.