

Product Name: Recombinant Human UROD (N-6His)
Catalog #: PEH1795




Summary

Name	UROD
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Uroporphyrinogen Decarboxylase is produced by our E.coli expression system and the target gene encoding Met1-Asn367 is expressed with a 6His tag at the N-terminus.
Accession #	AAH01778.1
Host	E.coli
Species	Human
Predicted Molecular Mass	43 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 1mM EDTA, pH 8.0.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image

Product Name: Recombinant Human UROD (N-6His)
Catalog #: PEH1795



Alternative Names

Uroporphyrinogen Decarboxylase; UPD; URO-D; UROD

Background

Uroporphyrinogen decarboxylase (UROD), is an enzyme of the heme biosynthetic pathway which belongs to the uroporphyrinogen decarboxylase family. This enzyme is responsible for catalyzing the conversion of uroporphyrinogen to coproporphyrinogen through the removal of four carboxymethyl side chains. UROD is a homodimeric enzyme that catalyzes the fifth step in heme biosynthesis: the elimination of carboxyl groups from the four acetate side chains of uroporphyrinogen III to yield coproporphyrinogen III. Defects in UROD are the cause of familial porphyria cutanea tarda (FPCT) and hepatoerythropoietic porphyria (HEP).

Note

For Research Use Only , Not for Diagnostic Use.