


Product Name: Recombinant Human UROIIIIS (C-6His)
Catalog #: PEH1797



Summary

Name	Uroporphyrinogen-III synthase/UROS
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Uroporphyrinogen-III Synthase is produced by our E.coli expression system and the target gene encoding Met1-Cys265 is expressed with a 6His tag at the C-terminus.
Accession #	P10746
Host	E.coli
Species	Human
Predicted Molecular Mass	29.7 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, pH 8.0.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	Uroporphyrinogen-III Synthase; UROIIIIS; UROS; Hydroxymethylbilane Hydrolyase [Cyclizing]; Uroporphyrinogen-III Cosynthase; UROS
Background	Uroporphyrinogen-III Synthase is an enzyme which belongs to the

Product Name: Recombinant Human UROIIIIS (C-6His)
Catalog #: PEH1797



uroporphyrinogen-III synthase family. Uroporphyrinogen-III Synthase is ubiquitous and it is involved in Porphyrin metabolism. Porphyrins act as cofactors for a multitude of enzymes that perform a variety of processes within the cell such as Methionine synthesis (Vitamin B12) or oxygen transport (Heme). Uroporphyrinogen-III Synthase can catalyze cyclization of the linear Tetrapyrrole, Hydroxymethylbilane, to the Macrocyclic Uroporphyrinogen III, the branch point for the various sub-pathways leading to the wide diversity of Porphyrins. Defects in Uroporphyrinogen-III Synthase are the cause of Congenital Erythropoietic Porphyria (CEP).

Note

For Research Use Only , Not for Diagnostic Use.