

Product Name: Recombinant Human VSNL1 (N-6His)
Catalog #: PEH1817




Summary

Name	Visinin-like protein 1/HLP3/VSNL1
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Visinin-Like Protein 1 is produced by our E.coli expression system and the target gene encoding Met1-Lys191 is expressed with a 6His tag at the N-terminus.
Accession #	P62760
Host	E.coli
Species	Human
Predicted Molecular Mass	24.3 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 20mM NaCl, pH 8.0.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human VSNL1 (N-6His)
Catalog #: PEH1817



Alternative Names

Visinin-Like Protein 1; VILIP; VLP-1; Hippocalcin-Like Protein 3; HLP3; VSNL1; VISL1

Background

Visinin-Like Protein 1 (VILIP) is a member of the Visinin/Recoverin subfamily of neuronal calcium sensor proteins. VILIP is strongly expressed in the Granule Cells of the Cerebellum where it associates with membranes in a Calcium-dependent manner and modulates intracellular signaling pathways of the central nervous system by directly or indirectly regulating the activity of Adenylyl Cyclase. It has been shown that VILIP regulates the inhibition of rhodopsin phosphorylation in a Calcium-dependent manner in vitro.

Note

For Research Use Only , Not for Diagnostic Use.