

Product Name: Recombinant Human WARS (N-6His)
Catalog #: PEH1829




Summary

Name	WARS/TrpRS
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Tryptophan-tRNA ligase, Cytoplasmic is produced by our E.coli expression system and the target gene encoding Met1-Gln471 is expressed with a 6His tag at the N-terminus.
Accession #	P23381
Host	E.coli
Species	Human
Predicted Molecular Mass	55.3 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 8% Trehalose, 4% Mannitol, 50mM NaCl, 0.05% Tween 80, pH7.5.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human WARS (N-6His)
Catalog #: PEH1829



Alternative Names

WARS also known as Tryptophanyl-tRNA synthetase; Interferon-induced protein 53

Background

There exists two types of tryptophanyl tRNA synthetases, the cytoplasmic form called WARS, the mitochondrial form called WARS2. WARS catalyzes the aminoacylation of tRNA (trp) with tryptophan and is induced by interferon. WARS regulates ERK, Akt, eNOS activation pathway, which are related with angiogenesis, cytoskeletal reorganization and shear stress-responsive gene expression.

Note

For Research Use Only , Not for Diagnostic Use.