

**Product Name: Recombinant Human Claudin-18.2 (N-8His)**  
**Catalog #: PEH2010**




---

## Summary

<b>Name</b>	CLDN18.2/Claudin 18.2
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	<1 EU/μg as determined by LAL test.
<b>Construction</b>	Recombinant Human CLDN18.2 is produced by our E.coli expression system and the target gene encoding Asp28-Leu76 is expressed with a 8His tag at the N-terminus.
<b>Accession #</b>	P56856-2
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Predicted Molecular Mass</b>	19.7 KDa
<b>Formulation</b>	Supplied as a 0.2 μm filtered solution of 25mM Tris-HCl, 25mM NaCl, 0.1% Triton X-100, 10% glycerol, 0.1mM GSH, 0.01mM GSSG, pH 8.0.
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Reconstitution</b>	

## SDS-PAGE image

Product Name: Recombinant Human Claudin-18.2 (N-8His)  
Catalog #: PEH2010



### Alternative Names

Claudin-18; CLDN18

### Background

Claudin-18 (CLDN18) is a protein that in humans is encoded by the CLDN18 gene. It belongs to the group of claudins. CLDN18 belongs to the large claudin family of proteins, which form tight junction strands in epithelial cells. CLDN18 plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity. CLDN18 has two isoform A1 and isoform A2. IMAB362 (Claudiximab) is a monoclonal antibody against isoform 2 of Claudin-18. It is under investigation for the treatment of gastrointestinal adenocarcinomas and pancreatic tumors. IMAB362 was developed by Ganymed Pharmaceuticals AG.

### Note

For Research Use Only , Not for Diagnostic Use.