



Product Name: Recombinant Human BIRC5 Protein
Catalog #: PEH2718



Summary

| | |
|------------------------------|--|
| Name | BIRC5_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-mSUMO & C-His |
| Accession # | O15392 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 24 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | AGFIHCPTEN EPDLAQCFEC FKELEGWEPD DDPIEEHKKH SSGCAFLSVK KQFEELTLGE FLKLDREKAK NKIAKETNNK KKEFEETAKK |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Baculoviral IAP repeat-containing protein 5 (Apoptosis inhibitor 4) (Apoptosis inhibitor survivin) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use